




TP53 Antibody

Product Code	CSB-RA825742A0HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04637
Immunogen	A synthesized peptide derived from human Mutant p53
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:50-1:200, IF:1:20-1:200
Relevance	<p>Acts as a tumor suppressor in many tumor types; induces growth arrest or apoptosis depending on the physiological circumstances and cell type. Involved in cell cycle regulation as a trans-activator that acts to negatively regulate cell division by controlling a set of genes required for this process. One of the activated genes is an inhibitor of cyclin-dependent kinases. Apoptosis induction seems to be mediated either by stimulation of BAX and FAS antigen expression, or by repression of Bcl-2 expression. In cooperation with mitochondrial PPIF is involved in activating oxidative stress-induced necrosis; the function is largely independent of transcription. Induces the transcription of long intergenic non-coding RNA p21 (lincRNA-p21) and lincRNA-Mkln1. LincRNA-p21 participates in TP53-dependent transcriptional repression leading to apoptosis and seem to have to effect on cell-cycle regulation. Implicated in Notch signaling cross-over. Prevents CDK7 kinase activity when associated to CAK complex in response to DNA damage, thus stopping cell cycle progression. Isoform 2 enhances the transactivation activity of isoform 1 from some but not all TP53-inducible promoters. Isoform 4 suppresses transactivation activity and impairs growth suppression mediated by isoform 1. Isoform 7 inhibits isoform 1-mediated apoptosis. Regulates the circadian clock by repressing CLOCK-ARNTL/BMAL1-mediated transcriptional activation of PER2 (PubMed:24051492).</p>
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Affinity-chromatography
Isotype	Rabbit IgG
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling; Cancer; Cell biology
Gene Names	TP53
Accession NO.	5C3


Image

Western Blot

Positive WB detected in: U-251 whole cell lysate

All lanes: Mutant p53 antibody at 1:1000


Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 38, 39, 40, 34, 35, 30, 24,

25 kDa

Observed band size: 53 kDa



IHC image of CSB-RA825742A0HU diluted at

1:100 and staining in paraffin-embedded human

glioma cancer performed on a Leica Bond™

system. After dewaxing and hydration, antigen

retrieval was mediated by high pressure in a

citrate buffer (pH 6.0). Section was blocked with


10% normal goat serum 30min at RT. Then

primary antibody (1% BSA) was incubated at

4°C overnight. The primary is detected by a Goat

anti-rabbit IgG polymer labeled by HRP and

visualized using 0.05% DAB.



Immunofluorescence staining of HepG2 Cells

with CSB-RA825742A0HU at 1:50, counter-

stained with DAPI. The cells were fixed in 4%

formaldehyde, permeated by 0.2% TritonX-100,

and blocked in 10% normal Goat Serum. The

cells were then incubated with the antibody

overnight at 4°C. Nuclear DNA was labeled in

blue with DAPI. The secondary antibody was

FITC-conjugated AffiniPure Goat Anti-Rabbit IgG

(H+L).