




Recombinant Human Zinc finger protein 384 (ZNF384), partial

Product Code	CSB-RP034554h
Relevance	Transcription factor that binds the consensus DNA sequence [GC]AAAAA. Seems to bind and regulate the promoters of MMP1, MMP3, MMP7 and COL1A1 (By similarity).
Abbreviation	Recombinant Human ZNF384 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TF68
Alias	CAG repeat protein 1 CAS-interacting zinc finger protein Nuclear matrix transcription factor 4 Short name: Nuclear matrix protein 4 Trinucleotide repeat-containing gene 1 protein
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MEESHFNSNPYFWPSIPTVSGQIENTMFINKMKDQLLPEKGCGLAPPHYPTLL TVPASVSLPSGISMDTESKSDQLTPHSQASVTQNITVVPVPSTGLMTAGVSCS QRWRREGSQSRGPGLVITSPSGSLVTTASSAQTFPISAPMIVSALPPGSQALQ VVPDLSKKVASTLTEEGGGGGGGGGSVAPKPPRGRKKKRMLESGLPEMNDP YVLSPEDDDDHQKDGTKYRCRMCSLTFYSKSEMQUI
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	ZNF384
Expression Region	1-245aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	53.1kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.