




Recombinant Human Protein phosphatase 1 regulatory subunit 11 (PPP1R11)

Product Code	CSB-RP038344h
Relevance	Inhibitor of protein phosphatase 1.
Abbreviation	Recombinant Human PPP1R11 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O60927
Alias	Hemochromatosis candidate gene V protein ;HCG VProtein phosphatase inhibitor 3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MAEAGAGLSETVTETT VTTTEPENRSLTIKLRKRKPEKKVEWTS DTV DNEHM GRRSSKCCCIYEKPRAFGESSTESDEEEEEEGCGHTHCVRGHRKGRRRATLG PTPTTPPQPPDPSQPPPGPMQH
Research Area	Others
Source	E.coli
Target Names	PPP1R11
Expression Region	1-126aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	41.0kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.