




Recombinant Human T-cell surface glycoprotein CD3 epsilon chain (CD3E), partial,Yeast

Product Code	CSB-YP004931HU
Relevance	The CD3 complex mediates signal transduction, resulting in T-cell activation and proliferation. Required for normal immune responses (PubMed:15546002, PubMed:8490660).
Abbreviation	Recombinant Human CD3E protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07766
Alias	T-cell surface antigen T3/Leu-4 epsilon chain;CD_antigen: CD3e
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Research Area	Immunology
Source	Yeast
Target Names	CD3E
Protein Names	Recommended name: T-cell surface glycoprotein CD3 epsilon chain Alternative name(s): T-cell surface antigen T3/Leu-4 epsilon chain CD_antigen= CD3e
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged Tag-Free
Mol. Weight	13.8kDa
Protein Length	Extracellular Domain

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.