



# Recombinant Human T-cell surface glycoprotein CD8 alpha chain (CD8A), partial

<b>Product Code</b>	CSB-YP004966HU
<b>Relevance</b>	Identifies cytotoxic/suppressor T-cells that interact with MHC class I bearing targets. CD8 is thought to play a role in the process of T-cell mediated killing. CD8 alpha chains binds to class I MHC molecules alpha-3 domains.
<b>Abbreviation</b>	Recombinant Human CD8A protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P01732
<b>Alias</b>	T-lymphocyte differentiation antigen T8/Leu-2 CD_antigen: CD8a
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	SQFRVSPLDRTWNLGETVELKCQVLLSNPTSGCSWLFQPRGAAASPTFLLYL SQNKPKAAEGLDTQRFSGKRLGDTFVLTLSDFRRENEGYYFCSALSNSIMYFS HFVPVFLPAKPTTTPAPRPPTPAPTIASQPLSLRPEACRPAAGGAVHTRGLDFA CD
<b>Research Area</b>	Immunology
<b>Source</b>	Yeast
<b>Target Names</b>	CD8A
<b>Protein Names</b>	Recommended name: T-cell surface glycoprotein CD8 alpha chain Alternative name(s): T-lymphocyte differentiation antigen T8/Leu-2 CD_antigen= CD8a
<b>Expression Region</b>	22-182aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	19.6kDa
<b>Protein Length</b>	Extracellular Domain
<b>Image</b>	

