




Recombinant Human Neuronal acetylcholine receptor subunit alpha-3 (CHRNA3), partial

Product Code	CSB-YP005389HU1
Abbreviation	Recombinant Human CHRNA3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32297
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	SEAEHRLFERLFEDYNEIIRPVANVSDPVIIHFEVSMSQLVKVDEVNQIMETNLW LKQIWNDYKLKWNPSDYGGAEFMRVPAQKIWKPDIVLYNNAVGFQVDDKTK ALLKYTG EVTWIPPAIFKSSCKIDVTYFPFDYQNCTMKFGSWSYDKAKIDLVLIG SSMNLKDYWESGEWAIKAPGYKHDIKYNCC EEIYPDITYSLYIRRL
Research Area	Neuroscience
Source	Yeast
Target Names	CHRNA3
Protein Names	Recommended name: Neuronal acetylcholine receptor subunit alpha-3
Expression Region	32-240aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	26.0 kDa
Protein Length	Extracellular Domain
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.