




Recombinant Human Collagen alpha-1 (XI) chain (COL11A1), partial

Product Code	CSB-YP005716HU
Relevance	May play an important role in fibrillogenesis by controlling lateral growth of collagen II fibrils.
Abbreviation	Recombinant Human COL11A1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P12107
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	GPMGLTGRPGPVGGPGSSGAKGESGDPGPQGPRGVQGPPGPTGKPGKRG RPGADGGRGMPGEPGAKGDRGFDGLPGLPGDKGHRGERGPQGPPGPPGD DGMRGEDGEIGPRGLPGEAGPRGLLGPRGTPGAPGQPGMAGVDGPPGPKG NMGPQGEPGPPGQQGNPGPQ
Source	Yeast
Target Names	COL11A1
Protein Names	Recommended name: Collagen alpha-1(XI) chain
Expression Region	532-699aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	42.8kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.