



Recombinant Human Glutamate receptor 3 (GRIA3), partial


Product Code	CSB-YP009900HU
Relevance	Receptor for glutamate that functions as ligand-gated ion channel in the central nervous syst and plays an important role in excitatory synaptic transmission. L- glutamate acts as an excitatory neurotransmitter at many synapses in the central nervous syst. Binding of the excitatory neurotransmitter L-glutamate induces a conformation change, leading to the opening of the cation channel, and thereby converts the chical signal to an electrical impulse. The receptor then desensitizes rapidly and enters a transient inactive state, characterized by the presence of bound agonist. In the presence of CACNG4 or CACNG7 or CACNG8, shows resensitization which is characterized by a delayed accumulation of current flux upon continued application of glutamate.
Abbreviation	Recombinant Human GRIA3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P42263
Alias	AMPA-selective glutamate receptor 3GluR-CGluR-K3Glutamate receptor ionotropic, AMPA 3 ;GluA3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	SLLGHYKWEKVFVLYDTERGFSILQAIMEAAVQNNWQVTARSVGNIKDVQEFR RIIEEMDRRQEKRYLIDCEVERINTILEQVVILGKHSRGYHYMLANL
Research Area	Neuroscience
Source	Yeast
Target Names	GRIA3
Protein Names	Recommended name: Glutamate receptor 3 Short name= GluR-3 Alternative name(s): AMPA-selective glutamate receptor 3 GluR-C GluR-K3 Glutamate receptor ionotropic, AMPA 3 Short name= GluA3
Expression Region	151-250aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	14.0kDa



Protein Length

Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.