




Recombinant Rabbit Interleukin 17A (IL17A), partial

Product Code	CSB-YP011597RB
Abbreviation	Recombinant Rabbit IL17A protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	G1SLF2
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	VKNGIAMPRNPGCPNAEDKNFPQNVKVSILNLSVNSRRPSDYNNRSTSPW TLHRNEDRERYPSVIWEAKCRHLGCVNAEGNEDHHMNSVPIQQEILVLRRES QHCPHSFRLEKMLVAVGCTCVTPIIHHMA
Source	Yeast
Target Names	IL17A
Expression Region	21-153aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	17.3 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.