




Recombinant Rabbit Interleukin-2 (IL2)

Product Code	CSB-YP011629RB
Abbreviation	Recombinant Rabbit IL2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O77620
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	APTSSSTKETQEQLDQLLLDLQVLLKGVNDYKNSKLSRMLTFKFKYMPKKVTEL KHLQCLEEEELKPLEEVLNLAQGKNSHGGNTRESISNINVTVLKLGSETFMCEY DETVTIVEFLNRWITFCQSIISASSS
Research Area	Cancer
Source	Yeast
Target Names	IL2
Protein Names	Recommended name: Interleukin-2 Short name= IL-2 Alternative name(s): T-cell growth factor Short name= TCGF
Expression Region	21-153aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	17.1kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.