




Recombinant Human Mucin-2 (MUC2), partial

Product Code	CSB-YP015222HU
Relevance	Coats the epithelia of the intestines, airways, and other mucus mbrane-containing organs. Thought to provide a protective, lubricating barrier against particles and infectious agents at mucosal surfaces. Major constituent of both the inner and outer mucus layers of the colon and may play a role in excluding bacteria from the inner mucus layer.
Abbreviation	Recombinant Human MUC2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q02817
Alias	Intestinal mucin-2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	VCSTWGNFHYKTFDGDVFRFPGLCDYNFASDCRGSYKEFAVHLKRGPGQAE APAGVESILLTIKDDTIYLTRHLAVLNGAVVSTPHYSPGLLIEKSDAYTKVYSRA GLTLMWNREDALMLELDTKFRNHTCGLCGDYNGLQSYSEFLSDGVLFSPLEF GNMQKINQPDVVCEDPEEEVAPASCSEHRAECERLLTAEAFADCQDL
Research Area	Signal Transduction
Source	Yeast
Target Names	MUC2
Protein Names	Recommended name: Mucin-2 Short name= MUC-2Alternative name(s): Intestinal mucin-2
Expression Region	36-240aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	24.8kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.