




Recombinant Mouse Protein disulfide-isomerase A2 (Pdia2)

Product Code	CSB-YP017719MOc7
Abbreviation	Recombinant Mouse Pdia2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	D3Z6P0
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Mouse Protein disulfide-isomerase A2(Pdia2)
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	<p>QGEPPGGPSEVLPEEPTGEEVPKEDGILVLNHRTLSLALQEHSALMVEFYAPW CGHCKELAPEYSKAAALLAAESAVVTLAKVDGPAEPELTKEFEVVGYP TLKFF QNGNRTNP E EYAGPKTAEGIAEWLRRRVGPSATHLEDEEGVQALMAKWDMV VIGFFQDLQGKDMATFLALAKDALDMTFGFTDQPQLFEKFGLTKD TVVLFKKF DEGRADFPVDKETGLDLGDL SRFLVIHSMHLVTEFNSQTS PKIFA AKILNHLLLF VNQTLAQHRELLTDFREAA PFRGQVLFVMVDVAADNSHVLNYFGLK AEEAP TLRLINVETTKKYAPTGVIAITAASVA AFCQAVLHGEIKHYLLS QEIPPDWDQGP VKTLVSKNF EQVAFDETKNVFVKFYAPWCSHCKEMAPAW EALAEKYKDREDI VIAELDATANELEAFSVLGYPTLKFFPAGPDRKVIDYKSTRDLETFSKFLDSGG HLPKEEPKEPAASAPEAQANSTLGPKEEL</p>
Research Area	Others
Source	Yeast
Target Names	Pdia2
Expression Region	21-527aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	57.8 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.