




Recombinant Human Pregnancy zone protein (PZP), partial

Product Code	CSB-YP019131HU
Relevance	Is able to inhibit all four classes of proteinases by a unique 'trapping' mechanism. This protein has a peptide stretch, called the 'bait region' which contains specific cleavage sites for different proteinases. When a proteinase cleaves the bait region, a conformational change is induced in the protein which traps the proteinase. The entrapped enzyme rains active against low molecular weight substrates (activity against high molecular weight substrates is greatly reduced). Following cleavage in the bait region a thioester bond is hydrolyzed and mediates the covalent binding of the protein to the proteinase.
Abbreviation	Recombinant Human PZP protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20742
Alias	C3 and PZP-like alpha-2-macroglobulin domain-containing protein 6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MKPEAELSVSSVYNLLTVKDLTNFPDNDVQEQEEQGHCPRPFFIHNGAIYVPL SSNEADIYSFLKGMGLKVFTNSKIRKPKSCSVIPSVSAGAVGQGYGAGLGVV ERPYVPQLGTYNVIPLNNEQSSGPVPETVRSYFPETWIWELVAVNSSGVAEVG VTVPDTITEWKAGAFCLSEDAGLGISSTASLRAFQPFFVELTMPYSVIRGEVFT LKATVLNLYLPKCIRAEGIEQKTFSSMTCASGANVSEQLSLKLPSNVVKESARA SFSVLGDILGSAMQNIQNLLQMPYGCGEQNMVLFAPNIYVLNLYLNETQQLTQEI KAKAVGYLITGYQRQLNYKHQDGSYSTFG
Research Area	Neuroscience
Source	Yeast
Target Names	PZP
Protein Names	Recommended name: Pregnancy zone protein Alternative name(s): C3 and PZP-like alpha-2-macroglobulin domain-containing protein 6
Expression Region	472-821aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	65.3kDa


Protein Length

Partial of Isoform 2

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.