




# Recombinant Mouse Tryptophan 2,3-dioxygenase (Tdo2)

<b>Product Code</b>	CSB-YP023351MO
<b>Relevance</b>	Incorporates oxygen into the indole moiety of tryptophan. Has a broad specificity towards tryptamine and derivatives including D- and L-tryptophan, 5-hydroxytryptophan and serotonin .
<b>Abbreviation</b>	Recombinant Mouse Tdo2 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P48776
<b>Alias</b>	Tryptamin 2,3-dioxygenase
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	MSGCPFAGNSVG YTLKNVSMEDNEEDRAQTGVNRASKGGLIYGNYLQLEKIL NAQELQSEVKGNKIHDEHLFIITHQAYELWFKQILWELDSVREIFQNGHVNDER NMLKVIARMHRVVVIFKLLVQQFSVLETMTALDFNDFREYLSPASGFQSLQFRL LENKIGVLQSLRVPYNRKHYRDNFGGDYNELLLKSESEQTLLQLVEAWLERTP GLEPNGFNFWGKFEKNILKGLEEEFLRIQAKTDSEEKEEQMAEFRKQKEVLLC LFDEKRHDYLLSKGERRLSYRALQGALMIYFYREEPRFQVPFQLTSLMDIDTL MTKWRYNHVCMVHRMLGTKAGTGGSSGYHYLRSTVSDRYKVFVDLNLSTY LVPRHWVPKMNPPIHKFLYTAEYSDSSYFSSDESD
<b>Research Area</b>	Others
<b>Source</b>	Yeast
<b>Target Names</b>	Tdo2
<b>Protein Names</b>	Recommended name: Tryptophan 2,3-dioxygenase Short name= TDO EC= 1.13.11.11 Alternative name(s): Tryptamin 2,3-dioxygenase Tryptophan oxygenase Short name= TO Short name= TRPO Tryptophan pyrrolase Tryptophan
<b>Expression Region</b>	1-406aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	49.8kDa
<b>Protein Length</b>	Full Length



## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.