




# Recombinant Sheep Protransforming growth factor alpha (TGFA), partial

|                          |                                                                                                                                                                                                                                                                                 |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP023445SH                                                                                                                                                                                                                                                                  |
| <b>Relevance</b>         | TGF alpha is a mitogenic polypeptide that is able to bind to the EGF receptor/EGFR and to act synergistically with TGF beta to promote anchorage-independent cell proliferation in soft agar                                                                                    |
| <b>Abbreviation</b>      | Recombinant Sheep TGFA protein, partial                                                                                                                                                                                                                                         |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| <b>Uniprot No.</b>       | P98135                                                                                                                                                                                                                                                                          |
| <b>Alias</b>             | EGF-like TGF Short name: ETGF TGF type 1                                                                                                                                                                                                                                        |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                             |
| <b>Immunogen Species</b> | Ovis aries (Sheep)                                                                                                                                                                                                                                                              |
| <b>Purity</b>            | Greater than 90% as determined by SDS-PAGE.                                                                                                                                                                                                                                     |
| <b>Sequence</b>          | ENSTSALSDPPVAAAVVSHFNDPCDSHTQFCFHGTCRFLQEEKPACVCHSG<br>YVGARCEHADLLAVVAASQKKQ                                                                                                                                                                                                   |
| <b>Research Area</b>     | Others                                                                                                                                                                                                                                                                          |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                           |
| <b>Target Names</b>      | TGFA                                                                                                                                                                                                                                                                            |
| <b>Protein Names</b>     | Recommended name: Protransforming growth factor alpha Cleaved into the following chain: 1. Transforming growth factor alpha Short name= 2. TGF-alpha<br>Alternative name(s): EGF-like TGF Short name= ETGF TGF type 1                                                           |
| <b>Expression Region</b> | 24-97aa                                                                                                                                                                                                                                                                         |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                             |
| <b>Tag Info</b>          | N-terminal GST-tagged                                                                                                                                                                                                                                                           |
| <b>Mol. Weight</b>       | 34.9kDa                                                                                                                                                                                                                                                                         |
| <b>Protein Length</b>    | Extracellular Domain                                                                                                                                                                                                                                                            |
| <b>Image</b>             |                                                                                                                                                                                                                                                                                 |



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.