




# Recombinant *Macaca mulatta* Uncharacterized protein (OSM)

<b>Product Code</b>	CSB-YP2949MOWf3
<b>Abbreviation</b>	Recombinant Rhesus macaque OSM protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	F7GF43
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	<i>Macaca mulatta</i> (Rhesus macaque)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	MASMAAMGSCSKEYRMLLGQLQKQTDLMQDTSRLDPYIRIQGLDIPKLREHC RESPGAFPSEETLRGLGRRGFLQTLNATLGRVLHRLADLEQLPKAQDLERS GLNIEDLEKLQMARPNVLGLRNNVYCMAQLLDNSDMTEPTKAGRGTPQPPTP TPTSDVFQRKLEGCSFLRGYHRFMHSVGRVFSKWGESPNRSRRHSPHQALR KGVRRTRPSRKGNRLMPRGQLPR
<b>Research Area</b>	Immunology
<b>Source</b>	Yeast
<b>Target Names</b>	OSM
<b>Expression Region</b>	1-231aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	C-terminal 6xHis-Myc-tagged
<b>Mol. Weight</b>	27.1
<b>Protein Length</b>	Full Length

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



## **Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.