





Recombinant Cavia porcellus Lipocalin Cav p 2.0101 (Lcncavp2)

Product Code	CSB-YP306654GU
Abbreviation	Recombinant Guinea pig Lcncavp2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P83508
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Cavia porcellus Lipocalin Cav p 2.0101 (Lcncavp2)
Immunogen Species	Cavia porcellus (Guinea pig)
Sensitivity	Not Test
Purity	Greater than 95% as determined by SDS-PAGE. Greater than 95% as determined by SEC-HPLC.
Sequence	DSIDYSKVPGNWRTIAIAADHVEKIEVNGELRAYFRQVDCTEGCDKISITFYTNT DGVCTEHTVVGARNGENDVYTVDYAGENTFQILCNSDDAFVIGSVNTDQNGQ TTKEVAIAAKRNFLTPEQEQKFQKAVQNAGIPLNIRYVIETDTCPD
Research Area	Others
Source	Yeast
Target Names	Lcncavp2
Expression Region	17-170aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	18.9 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



The purity of Lcncavp2 was greater than 95% as determined by SEC-HPLC

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.