





Recombinant Schistosoma mansoni Major egg antigen


Product Code	CSB-YP319453SXQ
Abbreviation	Recombinant Schistosoma mansoni Major egg antigen protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P12812
Alias	p40
Product Type	Recombinant Protein
Immunogen Species	Schistosoma mansoni (Blood fluke)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MSGGKQHNAVSIPVNREQRSFEKQRRDLLTGLEHGGGAHRGNSIAPYTEDW PSTVDNWIDSSWKRWDDDMRRLRRGMFALLPLDTFSIGILENPFALMHQMDR QIQDIRERMGSLDVPSTGSVNDFLKDAYEVEGEDGKVHFKVRFDAQGFAPQDIN VTSSENRVTVHAKKETTTDGRKCSREFCRMVQLPKSIDDSQLKCRMTDDGVL MLEAPVKVDQNQSLTLNESGQVAVRPKSDNQIKAVPASQALVAKGVHGLSYV DDGSGGKRLHVEVAVDPVYKPEDLFVNVDNRRVVVSGRHHKQKSDQHGRSS SFAEFSQSYAIPETVDPLSVSAQVVGNTLVLEAPLEKQHAITH
Research Area	Others
Source	Yeast
Protein Names	Recommended name: Major egg antigen Alternative name(s): p40
Expression Region	1-354aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	41.4kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of Yeast host and target protein, the LC-MS/MS Analysis result of CSB-YP319453SXQ could indicate that this peptide derived from Yeast-expressed *Schistosoma mansoni* (Blood fluke) N/A.



Based on the SEQUEST from database of Yeast host and target protein, the LC-MS/MS Analysis result of CSB-YP319453SXQ could indicate that this peptide derived from Yeast-expressed *Schistosoma mansoni* (Blood fluke) N/A.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.