




Recombinant Human cytomegalovirus

Uncharacterized protein UL128 (UL128)

Product Code	CSB-YP321481HWV
Abbreviation	Recombinant Human cytomegalovirus Uncharacterized protein UL128
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P16837
Product Type	Recombinant Protein
Immunogen Species	Human cytomegalovirus (strain AD169) (HHV-5) (HCMV)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MSPKDLTPFLTTLWLLLGHSRVPRVRAEECCFEINVNHPPERCYDFKMCNRFT VALRCPDGEVCYSPEKTAIEIRGIVTTMTHSLTRQVVHNKLTSCNYPYLYLEAD GRIRCGKVNDKAQYLLGAAGSVPRWINLEYDKITRIVGLDQYLESVKKHKRLD VCRAKMGYMLQ
Research Area	Microbiology
Source	Yeast
Target Names	UL128
Protein Names	Recommended name: Uncharacterized protein UL128
Expression Region	1-171aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	21.7kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.