





Recombinant Mouse Lymphocyte antigen 6G (Ly6g), partial

Product Code	CSB-YP333994MO
Abbreviation	Recombinant Mouse Ly6g protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P35461
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 90% as determined by SDS-PAGE. Greater than 90% as determined by SEC-HPLC.
Sequence	LECYNCIGVPPETSCNTTTCPFSDGFCVALEIEVIVDSHRSKVKSNLCLPICPTT LDNTEITGNAVNVKTYCCKEDLCNAAVPTGGSSWTMAG
Research Area	others
Source	Yeast
Target Names	Ly6g
Protein Names	Ly-6G.1
Expression Region	27-119aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	11.9kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



The purity of Ly6g was greater than 90% as determined by SEC-HPLC

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.