




# Recombinant Bovine Interferon alpha-G (IFNAG)

Product Code	CSB-YP344226BO
Relevance	Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.
Abbreviation	Recombinant Bovine IFNAG protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P49877
Alias	IFN-alpha7
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	CHLPHTHSLANRRVLTLLRQLRRVSPSSCLQDRNDFAFPQEALGGSQQLQKAQ AISVLHEVTQHTFQFFSVEGSAAVWDESLLDKLRDALDQQLTDLQFCLRQEEG LRGAPLLKEDSSLAVRKYFHRLTLYLQEKRHSPCAWEVVRAEVMRAFSSSTNL QERFRRKD
Research Area	Others
Source	Yeast
Target Names	IFNAG
Protein Names	Recommended name: Interferon alpha-G Alternative name(s): IFN-alpha7
Expression Region	24-189aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	21.2kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.