




Recombinant Mouse Major urinary protein 11 (Mup11), partial

Product Code	CSB-YP361419MO
Relevance	Binds pheromones that are released from drying urine of males. These pheromones affect the sexual behavior of males.
Abbreviation	Recombinant Mouse Mup11 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04938
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	EKINGEWHTIILASDKREKIEDNGNFRLFLEQIHVLENSLVLFHTVRDEECSEL SMVADKTEKAGEYSVTYDGFNTFTIPKTDYDNFLMAHLINEKDGETFQLMGLY GREPDLSSDIKERFAQLCEEHGILRENIIDLSNANRCLQARE
Research Area	Others
Source	Yeast
Target Names	Mup11
Protein Names	Recommended name: Major urinary protein 11
Expression Region	32-181 aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	19.6kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.