





Recombinant Mouse Serine protease inhibitor A3K (Serpina3k)

Product Code	CSB-YP362214MO
Abbreviation	Recombinant Mouse Serpina3k protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07759
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Mouse Serine protease inhibitor A3K(Serpina3k)
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 90% as determined by SDS-PAGE. Greater than 90% as determined by SEC-HPLC.
Sequence	FPDGTKEMDIVFHEHQDNGTQDDSLTLASVNTDFAFSLYKKLALKNPDTNIVFS PLSISAALALVSLGAKGKTMEIELEGLKFNLTTETPEADIHQGFGNLLQSLSQPED QDQINIGNAMFIEKDLQILAEFHEKTRALYQTEAFTADFQQPTEAKNLINDYVSN QTQGMIKELISELDERTLMVLVNYIYFKGKWKISFDPQDTFESEFYLDEKRSVK VPMMKMKLLTTRHFRDEELSCSVLELKYTGNASALLILPDQGRMQQVEASLQP ETLRKWRKTLFPSQIEELNLPKFSIASNYRLEEDVLPENGIGKEVFTEQADLSGIT ETKKLSVSQVVHKAFLDVAETGTEAAAATGVIGGIRKAILPAVHFNRPFVFIYH TSAQSILFMAKVNNPK
Research Area	Others
Source	Yeast
Target Names	Serpina3k
Expression Region	22-418aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	45.7 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



The purity of Serpina3k was greater than 90% as determined by SEC-HPLC

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.