




# Recombinant Mycobacterium tuberculosis Low molecular weight antigen MTB12 (mtb12)

<b>Product Code</b>	CSB-YP363739MVZ
<b>Relevance</b>	May play a role in the development of protective immune responses.
<b>Abbreviation</b>	Recombinant Mycobacterium tuberculosis mtb12 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P9WIN6
<b>Alias</b>	CFP-2 Low molecular weight protein antigen 2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	DPASAPDVPTAAQLTSLNLSLADPNVSFANKGSLVEGGIGGTEARIADHKLKKA AEHGDLPLSFSVTNIQPAAAGSATADVSVSGPKLSSPVTQNVTFVNQGGWML SRASAMELLQAAGN
<b>Research Area</b>	Others
<b>Source</b>	Yeast
<b>Target Names</b>	mtb12
<b>Protein Names</b>	Recommended name: Low molecular weight antigen MTB12 Alternative name(s): CFP-2 Low molecular weight protein antigen 2
<b>Expression Region</b>	49-168aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	14.1kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.