




Recombinant Legionella pneumophila Outer membrane protein MIP (mip),Yeast

Product Code	CSB-YP401660LLJ
Abbreviation	Recombinant Legionella pneumophila Outer membrane protein MIP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A5IGB8
Product Type	Recombinant Protein
Immunogen Species	Legionella pneumophila (strain Corby)
Purity	>85% (SDS-PAGE)
Sequence	ATDATSLATD KDKLSYSIGA DLGKNFKNQG IDVNPEAMAK GMQDAMSGAQ LALTEQQMKD VLNKFQKDLM AKRTAEFNKK ADENKVKGEA FLTENKNKPG VVVLPSGLQY KVINAGNGVK PGKSDTVTVE YTGRLIDGTV FDSTEKTGKP ATFQVSQVIP GWTEALQLMP AGSTWEIYVP SGLAYGPRSV GGPIGNETL IFKIHLISVK KSS
Source	Yeast
Target Names	mip
Protein Names	Recommended name: Outer membrane protein MIP EC= 5.2.1.8Alternative name(s): Macrophage infectivity potentiator Peptidyl-prolyl cis-trans isomerase Short name= PPIase Rotamase
Expression Region	21-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged Tag-Free
Protein Length	full length protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.