




Recombinant Feline herpesvirus 1 Envelope glycoprotein D (US6), partial

Product Code	CSB-YP5353FAJ
Abbreviation	Recombinant Feline herpesvirus 1 Envelope glycoprotein D, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q89634
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Immunogen Species	Feline herpesvirus 1 (FeHV-1) (Feline viral rhinotracheitis virus)
Sensitivity	Not Test
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	LTTTPKTTTVYVKGFNIPPLRYNYTQARIVPKIPQAMDPKITAEVRYVTSMDSCG MVALISEPDIDATIRTIQLSQKKTYNATISWFKVTQGCEYPMFLMDMRLCDPKR EFGICALRSPSYWLEPLTKYMFLTDDDELGLIMMAPAQFNQGQYRRVITIDGSMF YTDFMVQLSPTPCWFAPDRYEEILHEWCNRNVKTIGLDGARDYHYYWVPYNP QPHHKAVLLYWYRTHGREPPVRFQEAIKYDRPAIPSGSEDSKRSNDSRGESS GPNWIDIENYTPKNNVPIISDDDVPTAPPKGMNNQS
Research Area	Others
Source	Yeast
Target Names	US6
Expression Region	25-328aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	36.7 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.