




# Recombinant Bovine Pregnancy-associated glycoprotein 2 (PAG2)

<b>Product Code</b>	CSB-YP637747BO
<b>Abbreviation</b>	Recombinant Bovine PAG2 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q28057
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	KKMKTLETLREKNLLNLFLEEQAYRLSKNDSKITIHPLRNYLDTAYVGNITIGT PPQEFRVVFDTGSANLWVPCITCTSPACYTHKTFNPQNSSSFREVGGSPITIFYG SGIIQGFLGSDTVRIGNLVSPEQSFGLSLEEYGFDSLPGILGLAFPAMGIEDTI PIFDNLWSHGAFSEPVFAFYLNNTNKPEGSVVMFGGVDHRYKYGELNWIPVSQ TSHWQISMNNISMNGTVTACSCGCEALLDTGTSMIYGPTKLVTNIHKLMNARL ENSEYVVSCDAVKTLPPVIFNINIDYPLRPQAYIIKIQNSCRSVFQGGTENS NTWILGDIFLRQYFSVFDRKNRRIGLAPAV
<b>Research Area</b>	Signal Transduction
<b>Source</b>	Yeast
<b>Target Names</b>	PAG2
<b>Protein Names</b>	Recommended name: Pregnancy-associated glycoprotein 2 Short name= PAG 2 EC= 3.4.23.-
<b>Expression Region</b>	22-376aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	42.1
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.