




# Recombinant Mouse Uncharacterized aarF domain-containing protein kinase 2 (Adck2), partial, Yeast

|                          |                                                                                                                                                                                                                                                                                 |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP761222MO                                                                                                                                                                                                                                                                  |
| <b>Relevance</b>         | The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).                                                                                                           |
| <b>Abbreviation</b>      | Recombinant Mouse Adck2 protein, partial                                                                                                                                                                                                                                        |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| <b>Uniprot No.</b>       | Q6NSR3                                                                                                                                                                                                                                                                          |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                             |
| <b>Immunogen Species</b> | Mus musculus (Mouse)                                                                                                                                                                                                                                                            |
| <b>Purity</b>            | Greater than 90% as determined by SDS-PAGE.                                                                                                                                                                                                                                     |
| <b>Research Area</b>     | Others                                                                                                                                                                                                                                                                          |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                           |
| <b>Target Names</b>      | Adck2                                                                                                                                                                                                                                                                           |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                             |
| <b>Tag Info</b>          | N-terminal His-tagged<br>Tag-Free                                                                                                                                                                                                                                               |
| <b>Mol. Weight</b>       | 47.4kDa                                                                                                                                                                                                                                                                         |
| <b>Protein Length</b>    | Partial                                                                                                                                                                                                                                                                         |

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### **Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.