




Recombinant Mouse Protein cereblon (Crbn)

Product Code	CSB-YP806030MO
Relevance	Substrate recognition component of a DCX (DDB1-CUL4-X-box) E3 protein ligase complex that mediates the ubiquitination and subsequent proteasomal degradation of target proteins, such as MEIS2. Normal degradation of key regulatory proteins is required for normal limb outgrowth and expression of the fibroblast growth factor FGF8. May play a role in memory and learning by regulating the assembly and neuronal surface expression of large-conductance calcium-activated potassium channels in brain regions involved in memory and learning via its interaction with KCNT1.
Abbreviation	Recombinant Mouse Crbn protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8C7D2
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MAGEGDQQDAAHNMGNHLPLLPADSEDEDDEIEMEVEDQDSKEARKPNIINF DTSLPTSHTYLGADMEEFHGRTLHDDSDCQVIPVLPEVLMILIPGQTLPLQLSH PQEVSMVRNLIQKDRTFAYLAYSNNVQEREAQFGTTAEIYAYREEQEFGIEVVK VKAIGRQRFKVLELRTQSDGIQQAQVQILPECVLPSTMSAVQLESLNKCQVFPS KPISWEDQYSCKWWQKYQKRKFHCANLTSPRWLYSLYDAETLMDRIKKQL REWDENLKDDSLPENPIDFSYRVAACLPIDVLRILQLLKIGSAIQRRLRCELDINN KCTSLCCKQCQETEITTKNEIFSLSLCGPMAAYVNPBGYVHETLTVYKASNLNL IGRPSTVHSWFPGYAWTIAQCKICASHIGWKFTATKKDMSPQKFWGLTRALL PTIPETEDEISPDKVILCL
Research Area	Neuroscience
Source	Yeast
Target Names	Crbn
Protein Names	Recommended name: Protein cereblon Short name= Protein PiL
Expression Region	1-445aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	54.9kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.