




Recombinant Mouse Histone deacetylase complex subunit SAP130 (Sap130), partial

Product Code	CSB-YP807359MO
Relevance	Acts as a transcriptional repressor. May function in the assembly and/or enzymatic activity of the mSin3A corepressor complex or in mediating interactions between the complex and other regulatory complexes
Abbreviation	Recombinant Mouse Sap130 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8BIH0
Alias	130 kDa Sin3-associated polypeptide Sin3-associated polypeptide p130
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	PRKQQHVISTEEGDMMETNSTDDEKSAKSLLVKAEKRKSPPEYIDEEGVRY VPVRPRPPITLLRHYNPNWKAAYHHFQRYSDVRVKEEKKAMLQEIANQKGVS CRAQGKWVHLCAAQLQLTNLEHDVYERLTNLQEGIIIPKKKAATDDDLHRINEL IQGNMQRCKLVMDQISEARDSMLKVLDPKDRVLKLLNKNKGTGKKVSKLKRKEK V
Research Area	Others
Source	Yeast
Target Names	Sap130
Protein Names	Recommended name: Histone deacetylase complex subunit SAP130 Alternative name(s): 130 kDa Sin3-associated polypeptide Sin3-associated polypeptide p130
Expression Region	845-1057aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	26.8kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.