




Recombinant Human Ribonucleoprotein PTB-binding 2 (RAVER2), partial

Product Code	CSB-YP872547HU
Relevance	May bind single-stranded nucleic acids.Curated
Abbreviation	Recombinant Human RAVR2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HCJ3
Alias	Protein raver-2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MAAAAGDGGGEGGAGLGSAAGLGPGPLRGQGPSAEAHEGAPDMPAALH PEEVAARLQRMQRELSNRRKILVKNLPQDSNCQEVHDLKDYDLKYCYVDRN KRTAFVTLNGEQAQNAIQMFHQYSFRGKDLIVQLQPT
Research Area	Epigenetics and Nuclear Signaling
Source	Yeast
Target Names	RAVER2
Protein Names	Recommended name: Ribonucleoprotein PTB-binding 2 Alternative name(s): Protein raver-2
Expression Region	1-140aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	17.0kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.