




Recombinant Human Plexin-A1 (PLXNA1), partial

Product Code	CSB-YP883448HU
Relevance	Coreceptor for SA3A, SA3C, SA3F and SA6D. Necessary for signaling by class 3 saphorins and subsequent rodeling of the cytoskeleton. Plays a role in axon guidance, invasive growth and cell migration. Class 3 saphorins bind to a complex composed of a neuropilin and a plexin. The plexin modulates the affinity of the complex for specific saphorins, and its Cytoplasmic domain is required for the activation of down-stream signaling events in the cytoplasm .
Abbreviation	Recombinant Human PLXNA1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UIW2
Alias	Semaphorin receptor NOV
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	LNAGSDVAVSVGGRPCSFSWRNSREIRCLTPPGQSPGSAPIIININRAQLTNPE VKYNYTEDPTILRIDPEWSINSGGTLLTVTGTNLATVREPRIRAKYGGIERENG LVYNDTTMVCRAPSVANPVRSPPELGERPDELGFVMDNVRSLVLNSTSFLYY PDPVL
Research Area	Neuroscience
Source	Yeast
Target Names	PLXNA1
Protein Names	Recommended name: Plexin-A1 Alternative name(s): Semaphorin receptor NOV
Expression Region	986-1152aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	20.3kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.