



TUBE CHECKER (Black)

Catalog No. DAI-TC-BLACK

It is particularly useful for marking on plastic cassette of paraffin embedding.

TUBE CHECKER is a water/solvent proof twin-point pen.

- * "Tube Checker" pen enables you to write on paper, plastic, glass and metal.
- * The writings are waterproof after being dried and will not blur with water.
- * After drying, written parts become solvent proof against organic solvent mainly used in the lab such as ethanol, xylene and chloroform.

How to use:

- * Uncap the pen and start writing. The ink will dry in 20 seconds or so.
- * To avoid running, do not touch the ink with wet hands or dip in water before it is dried.

Color: Black

Caution:

- * Be sure to cap the pen after use.
- * Prevent direct exposure to the sun and keep it in a dark place.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN

[International/Outside Japan]

[国内連絡先]

TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618

TEL: 03-5632-9610 / FAX: 03-5632-9619

e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com

e-mail: support@cosmobio.co.jp / URL: www.cosmobio.co.jp