

Purified Antibody Data Sheet

Product name: Mouse anti-human CD63 antibody, 100µg

Product code : HAK-Anti-CD63-A1-9A2

Clone No. A1-9A2

Isotype mouse IgG1/kappa

Concentration 1mg/mL in PBS

Volume 100µL (100µg)

Storage Avoid freezing and thawing.

Store below -70°C upon receipt.

Small aliquots are recommended when using.

Immunogen Protein expressed by transfected plasmid DNA (CD63)

Purification Purified from hybridoma culture supernatant (serum-free medium) using a protein G column.

Lot No. Indicated on the tube label.

Application WB (see attachment), ELISA

Western Blot

SDS-PAGE was performed with 5 µg of cell extract of the immortalized human mesenchymal stem cell line UE6E7T-3, and Western blotting was performed using this antibody (mouse anti-human CD63 monoclonal antibody, clone # A1-9A2) as the primary antibody and HRP-labeled anti-mouse IgG antibody as the secondary antibody.

