

Specification Sheet

Monoclonal Antibody Against Human Glycophorin A

Code	2AHBGPA
Clone No.	Hyb- 5399
Manufacturing Process	Hybridoma obtained by fusing mouse myeloma cells (NS-1) with BALB/c mouse spleen cells immunized with human erythrocyte is proliferated in the BALB/c mouse abdominal cavity and ascites is collected.
Specificity	This monoclonal antibody specifically reacts with human glycophorin A. (Hemagglutinated with all of ABO blood group and MN blood group.)
Purification	Ammonium sulfate precipitation.
Ig Class	IgM
Form	Saline (containing 0.1% NaN ₃)
Protein Concentration	0.5 mg/mL
Package	0.2 mL
Storage	2-8°C (Below -20°C for long term storage)
Application	ELISA
Reference	1) 岡本宏明, 他: 医学と生物学 112 : 145-148, 1986.



INSTITUTE OF IMMUNOLOGY Co., LTD.

1-1-10, Koraku, Bunkyo-ku, Tokyo, JAPAN 112-0004

Tel : +81-3-3814-4081 Fax : +81-3-3814-5957

e-mail : info@tokumen.co.jp

<http://www.tokumen.co.jp/english/>

2AHBGPA1210-03

This product is for *in vitro* research use only and is not intended for use humans or animals.

2012.12.14