

Specification Sheet

Monoclonal Antibody Against Human HSP90 alpha beta

Code	2AHSP20
Clone No.	Hyb- K3720A
Manufacturing Process	Hybridoma obtained by fusing mouse myeloma cells (NS-1) with BALB/c mouse spleen cells immunized with recombinant human Heat-Shock Protein (HSP) 90 beta expressed in E.coli is proliferated in the BALB/c mouse abdominal cavity and ascites is collected.
Specificity	This monoclonal antibody specifically reacts with conserved region of human HSP90 alpha and beta (amino acids 291-304).
Purification	Ammonium sulfate precipitation.
Ig Class	IgG ₁
Form	Saline (containing 0.1% NaN ₃)
Protein Concentration	1 mg/mL
Package	0.2 mL
Storage	2-8°C (Below -20°C for long term storage)
Application	ELISA, Western blotting, IF
Reference	1) Matsusaka T: Jpn J Oral Biol 38 : 633-646, 1996. 2) Nemoto T, et al: J Biol Chem 272 26179-26187, 1997. 3) Nemoto T, et al: Biochem Mol Biol Int 42 : 881-889, 1997. 4) Kishimoto J, et al: Cell Stress & Chaperones 10 : 296-311, 2005.



INSTITUTE OF IMMUNOLOGY Co., Ltd.

1-1-10, Koraku, Bunkyo-ku, Tokyo, JAPAN 112-0004

Tel : +81-3-3814-4081 Fax : +81-3-3814-5957

e-mail : info@tokumen.co.jp

<http://www.tokumen.co.jp/english/>

2AHSP201210-05

This product is for *in vitro* research use only and is not intended for use humans or animals.

2012.12.14