

Specification Sheet

Monoclonal Antibody Against Mouse RBP-J kappa

Code	2ARBP3
Clone No.	Hyb- T6719
Manufacturing Process	Hybridoma obtained by fusing mouse myeloma cells (NS-1) with PVG/cSea rat spleen cells immunized with mouse RBP-J kappa expressed in E.coli is cultured and culture supernatant is collected.
Specificity	This monoclonal antibody specifically reacts with RBP-J kappa.
Purification	Centrifugation
Ig Class	IgG _{2a}
Form	Culture supernatant containing NaN ₃ (culture medium : ASF-104).
Protein Concentration	
Package	1 mL
Storage	2-8°C (Below -20°C for long term storage)
Application	ELISA, Western blotting ^{2),3)}
Reference	1) Sakai T, et al: J Biochem 118 : 621-628, 1995. 2) Hamaguchi Y, et al: J Biochem 112 : 314-320, 1992. 3) Grego-Bessa J, et al: Developmental Cell 12 : 415-429, 2007.



INSTITUTE OF IMMUNOLOGY Co., Ltd.

1-1-10, Koraku, Bunkyo-ku, Tokyo, JAPAN 112-0004

Tel : +81-3-3814-4081 Fax : +81-3-3814-5957

e-mail : info@tokumen.co.jp

<http://www.tokumen.co.jp/english/>

2ARBP31210-03

This product is for *in vitro* research use only and is not intended for use humans or animals.

2012.12.14