

Datasheet for PNA – Peanut Lectin

PNA – Peanut Lectin	
PNA	Catalog Number: MGC-J114 Pkg: 50mg Product Form: Lyophilized powder (salt-free) Storage: Stable for 1 year at 4°C

Origin Derived from red kidney bean (*Phaseolus vulgaris*)

Sugar Specificity D-GalNAc

Mitogenic Activity Present

Blood Type Specificity Non-specific

Features This lectin fraction is enriched in isolectins from red kidney bean containing a higher proportion of L subunits (E4, E3L1, E2L2, E1L3, L4), which are known for their lymphocyte mitogenic activity.