

Datasheet for Solvent Set (for ABEE Labeling Kit)

Catalog No.J712 | Package: 1 SET

Contents:

- Borate buffer powder (29g) × 3 bottles
- Trifluoroacetic acid solution (6 mL) × 1 bottle

(Sufficient for approximately 80 analyses)

Storage:Store at room temperature (25°C)

Feature:

Allows easy preparation of the following HPLC mobile phases:

- Solvent A: 0.2M Potassium borate buffer (pH 8.9)/acetonitrile = 93/7
- Solvent B: 0.02% Trifluoroacetic acid/acetonitrile = 50/50

Other Required Reagents:HPLC-grade acetonitrile

References:

1. Yasuno, S. et al., Biosci. Biotechnol. Biochem., 61, 1944 (1997)
2. Yasuno, S. et al., Biosci. Biotechnol. Biochem., 63, 1353 (1999)