

Perseus Proteomics Inc. 4-7-6, Komaba, Meguro-ku, Tokyo 153-0041 Japan TEL: +81-3-5738-1705 FAX: +81-3-3481-5760 https://www.ppmx.com info@ppmx.com

Anti human HNF4 alpha mouse monoclonal antibody(ProteinA purified)

Hepatocyte nuclear factor 4 alpha

		Application	Application / Recommended Concentration		
Code No	PP-H1415-0C	In order to obtain the best results, optimal working dilutions should be determined by each individual user.			
Clone No.	H1415	Western	n Blot	1 ug/mL	
Lot.	B-1	Non rec	lucing Western Blot	Decide by use	
Concentration	1 mg/mL				
Volume	100 uL	ELISA		0.3 ug/mL	
lg Class	G2a	Immund	oprecipitation	Decide by use	
Description	Hepatocyte nuclear factor 4 alpha (HNF4, HNF4a; NR2A1) is a member of orphan nuclear receptor. HNF4a is expressed in the liver, kidney, intestine and	Supersl	nift Assay	2ng/assay	
	pancreas. Mutation of HNF4a in humans has been associated with maturity-onset diabetes of the young type 1 (MODY1). HNF4 binds to DNA as an exclusive homodimer. The HNF4a gene is alternatively spliced	Chroma	atin immunoprecipitatic	Decide by use	
	and may generate up to nine different isoforms, HNF4a1 through HNF4a9.	Immund	bhistochemistry	10-20 ug/mL	
Nomenclature	NR2A1			· SSER	
Genbank	X87871		Human Liver Hepatocyte Paraffin section	Human intestinum tenue	
Origin	Produced in serum-free medium with hybridoma of mouse myeloma cells (NS-1) and spleen cells derived from a BALB/c mouse immunized with Baculovirus- expressed recombinant human HNF4 alpha (393-461 aa)	Storage	Store at 2 - 8 °C up to one month. For long-term storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in a frost-free freezer is not recommended.		
Specificity	This antibody specifically recognizes human HNF4 alpha 1, 2, 4, 5, 7,8 and cross reacts with mouse and rat HNF4 alpha 1, 2, 4, 5, 7,8.	Reference	Tanaka T, <i>et al.</i> J Pathol. 2006; 208(5): 662-72 Sugai M, <i>et al.</i> Pathol Int. 2008; 58(11): 681-6		
Purification	Affinity chromatography with Protein A				
Formulation	Physiological saline with 0.1% NaN3 as a preservative.	Notes	,	eact with lead and copper plumbing etal azides. Flush with large amounts sal.	

FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

Not for Diagnostic or Therapeutic use. Purchase of this product does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written consent of Perseus Proteomics Inc. is prohibited.

MADE IN JAPAN

Mar 9, 2022