


Storage box for Fecal Collection Samplers

Catalog numbers: TSG-FS-0011B, TSG-FS-0011M, TSG-FS-0011S

- This product is a paper storage box specifically designed for fecal collection samplers.
- All samplers can be stored upright, and the box is supplied with a lid.
- It is made of a material suitable for freezer use (down to -30°C), and essential information can be written directly on the surface.

Please note: Brush-type samplers should not be stored in the freezer, as this may damage the container.



Catalog Number	Product Name	Capacity	Dimensions
TSG-FS-0011B	Storage box (for brush type Feces Sampler)	100 Samplers	L147 mm x W147 mm x H125 mm
TSG- FS-0011M	Storage Box for Metabolo Keeper®	16 Samplers	L147 mm x W147 mm x H125 mm
TSG- FS-0011S	Storage box (for spoon type Feces Sampler)	25 Samplers	L147 mm x W147 mm x H97 mm