

# LC3b [6His-tagged]

## Modifying Protein

Alternate Names: Microtubule-Associated Protein 1, Light Chain 3, Beta; MAP1LC3B

Cat. No. 60-0110-500

Lot. No. 30097

Quantity: 500 µg

Storage: -70°C

FOR RESEARCH USE ONLY

NOT FOR USE IN HUMANS



CERTIFICATE OF ANALYSIS Page 1 of 2

## Background

The enzymes of the autophagy pathway play a pivotal role in the degradation of cytoplasmic constituents and organelles. Structures known as autophagosomes sequester portions of the cytoplasm which are degraded by the lysosome and recycled back into the cell. (Kuma *et al.* 2004). Three classes of enzymes are involved in autophagy; E1-like activating enzymes (E1s), E2-like conjugating enzymes (E2s) and ubiquitin-like proteins (ubls). Microtubule Associated Protein 1, Light Chain 3 beta, MAP1LC3b or LC3b is a member of the Ubl family and cloning of the human gene was first described by He *et al.* (2003). The C-terminal of LC3b is cleaved by ATG4B C-terminal to Gly120 *in vitro* and this has been shown to be required for the formation of intermediates with the ubiquitin-like activating enzyme ATG7 (Tanida *et al.* 2004). A proteomic analysis of the autophagosome interaction network was performed on a human cell line undergoing a basal level of autophagy. This study revealed that six of the ATG8 homologues (including LC3b) interacted with a cohort of 67 proteins in which there was frequent involvement with a conserved ATG8 surface region shown previously to interact with LC3 interacting regions in partner proteins (Behrends *et al.* 2010).

### References:

Behrends C, Sowa ME, Gygi SP, Harper JW (2010) Network organization of the human autophagy system. *Nature* 466, 68-76.

Continued on page 2

## Physical Characteristics

**Species:** human

**Source:** *E. coli*

**Quantity:** 500 µg

**Concentration:** 1 mg/ml

**Formulation:** 50 mM HEPES pH 7.5, 150 mM sodium chloride, 2 mM dithiothreitol, 10% glycerol

**Molecular Weight:** 16.93 kDa

**Purity:** >98% by InstantBlue™ SDS-PAGE

**Stability/Storage:** 12 months at -70°C; aliquot as required

### Protein Sequence:


**MGSSHHHHHSSGENLYFQGHMLGSMPSSEK**  
TFKQRRTFEQRVEDVRLIREQHPTKIPVI  
IERKGEKQLPVLDTKFLVDPDHVNMSELIKI  
IRRLQLNANQAFFLLVNGHSMVSVSTPI  
SEVYESEKDEDDGFLYVMVYASQETFG

Tag (**bold text**): N-terminal His  
Protease cleavage site: TEV (**ENLYF**▼**QG**)  
LC3b (regular text): Start **bold italics** (amino acid residues 1-120)  
Accession number: NP\_073729.1

## Quality Assurance

### Purity:

4-12% gradient SDS-PAGE  
InstantBlue™ staining  
Lane 1: MW markers  
Lane 2: 1 µg His-LC3b



### Protein Identification:

Confirmed by mass spectrometry.

### E1 Thioester LC3b Loading Assay:

The activity of LC3b was validated by loading LC3b onto the active cysteine of His-ATG7. Incubation of the His-ATG7 enzyme in the presence of LC3b and ATP at 30°C was compared at two time points, T<sub>0</sub> and T<sub>10</sub> minutes. Sensitivity of the LC3b/His-ATG7 thioester bond to the reducing agent DTT was confirmed.



www.ubiquigent.com  
Dundee, Scotland, UK

### ORDERS / SALES SUPPORT

International: +1-617-245-0003  
US Toll-Free: 1-888-4E1E2E3 (1-888-431-3233)  
Email: sales.support@ubiquigent.com

### UK HQ and TECHNICAL SUPPORT

International: +44 (0) 1382 381147 (9AM-5PM UTC)  
US/Canada: +1-617-245-0020 (9AM-5PM UTC)  
Email: tech.support@ubiquigent.com

Email services@ubiquigent.com for enquiries regarding compound profiling and/or custom assay development services.

© Ubiquigent 2012. Unless otherwise noted, Ubiquigent, Ubiquigent logo and all other trademarks are the property of Ubiquigent, Ltd.

Limited Terms of Use: For research use only. Not for use in humans or for diagnostics. Not for distribution or resale in any form, modification or derivative OR for use in providing services to a third party (e.g. screening or profiling) without the written permission of Ubiquigent, Ltd.

Lot-specific COA version tracker: v1.0.0

# LC3b [6His-tagged]

## Modifying Protein

**Alternate Names:** Microtubule-Associated Protein 1, Light Chain 3, Beta; MAP1LC3B

**Cat. No.** 60-0110-500

**Lot. No.** 30097

**Quantity:** 500 µg

**Storage:** -70°C

FOR RESEARCH USE ONLY

NOT FOR USE IN HUMANS



CERTIFICATE OF ANALYSIS Page 2 of 2

## Background

*Continued from page 1*

He H, Dang Y, Dai F, Guo Z, Wu J, et al. (2003) Post-translational modifications of three members of the human MAP1LC3 family and detection of a novel type of modification for MAP1LC3B. *J Biol Chem* **278**, 29278-29287.

Kuma A, Hatano M, Matsui M, Yamamoto A, Nakaya H, et al. (2004) The role of autophagy during the early neonatal starvation period. *Nature* **432**, 1032-1036.

Tanida I, Ueno T, Kominami E (2004) Human light chain 3/MAP1LC3B is cleaved at its carboxyl-terminal Met121 to expose Gly120 for lipidation and targeting to autophagosomal membranes. *J Biol Chem* **279**, 47704-47710.



[www.ubiquigent.com](http://www.ubiquigent.com)  
Dundee, Scotland, UK

### ORDERS / SALES SUPPORT

**International:** +1-617-245-0003  
**US Toll-Free:** 1-888-4E1E2E3 (1-888-431-3233)  
**Email:** [sales.support@ubiquigent.com](mailto:sales.support@ubiquigent.com)

### UK HQ and TECHNICAL SUPPORT

**International:** +44 (0) 1382 381147 (9AM-5PM UTC)  
**US/Canada:** +1-617-245-0020 (9AM-5PM UTC)  
**Email:** [tech.support@ubiquigent.com](mailto:tech.support@ubiquigent.com)

Email [services@ubiquigent.com](mailto:services@ubiquigent.com) for enquiries regarding compound profiling and/or custom assay development services.

© Ubiquigent 2012. Unless otherwise noted, Ubiquigent, Ubiquigent logo and all other trademarks are the property of Ubiquigent, Ltd.

**Limited Terms of Use:** For research use only. Not for use in humans or for diagnostics. Not for distribution or resale in any form, modification or derivative OR for use in providing services to a third party (e.g. screening or profiling) without the written permission of Ubiquigent, Ltd.

Lot-specific COA version tracker: v1.0.0