



YAMASA CORPORATION

DIAGNOSTICS DIVISION
Tokyo Office
1-23-8, Nihonbashi-kakigaracho
Chuo-ku, Tokyo, 103-0014, Japan
Tel. +81-3-3668-8558 Fax. +81-3-3668-8407
<https://diagnostics.yamasa.com/>

Anti-Human IgG3 Monoclonal Antibody

Abbreviation : MCA CG-01-DE

Code No. : 7645

Quantity : 100 μ L (100 μ g)

Clone : 5H9-2a

Quality : BALB/c IgG2b κ (Purified by DEAE chromatography.)

Solution : 1 mg/mL IgG, 0.1% NaN₃, PBS, pH 7.4

Character :	Reactivity		Reactivity	
to	Human IgA	< 0.5%	Human IgG2	< 0.5%
	Human IgE	< 0.5%	Human IgG4	< 0.5%
	Human IgG1	< 0.5%	Human IgM	< 0.5%

Storage: The solution should be stored below -20°C.

Repeated freezing and thawing should be avoided.

Note : Research use only.



COSMO BIO Co., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
e-mail: support@cosmobiousa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobiousa.com