



YAMASA CORPORATION

DIAGNOSTICS DIVISION
Tokyo Office
1-23-8, Nihonbashi-kakigaracho
Chuo-ku, Tokyo, 103-0014, Japan
Tel. +81-3-3668-8558 Fax. +81-3-3668-8407
<https://diagnostics.yamasa.com/>

Anti-Human Cardiac Myosin Light Chain I Monoclonal Antibody (520)

Code No. : 7887
Product name: ANTI-HUMAN CARDIAC MYOSIN LIGHT CHAIN I(520)
Quantity: 200 μ L (200 μ g)
Clone: MLM520
Quality: Balb/c IgG1 (Purified by protein A chromatography).
Solution: 1mg/mL IgG, 0.1%NaN₃, PBS, pH7.4
Character: This antibody reacts human, porcine, bovine, rat ventricle myosin light chain I.
It is possible to use for western blotting.
Storage: The solution should be stored below -20°C.
Repeated freezing and thawing should be avoided.
Note: Research use only.



COSMO BIO Co., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
e-mail: support@cosmobiousa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobiousa.com

Feb. 26, 2026