



YAMASA CORPORATION

DIAGNOSTICS DIVISION
Tokyo Office
1-23-8, Nihonbashi-kakigaracho
Chuo-ku, Tokyo, 103-0014, Japan
Tel. +81-3-3668-8558 Fax. +81-3-3668-8407
<https://diagnostics.yamasa.com/>

Anti-Mouse Interferon β Monoclonal Antibody

Code No. : 7891
Quantity: 500 μ L (1 mg)
Clone: 7F-D3
Quality: ACI rat IgG1 (Purified by Protein G chromatography.)
Solution: 2 mg/mL IgG, PBS, pH 7.4
Character:

Reactivity

to	Mouse IFN β	(L cell)	+
		(Recombinant)	+
		α , γ	—
	Human IFN		—

Example: Plate 3×10^4 cells/100 μ L of L cells in MEM containing 10% calf serum to a well of a micro plate.

Incubate at 37 °C for 24 hours.

Add 50 μ L of diluted 7891 and 50 μ L of mouse IFN- β in MEM.

Incubate at 37 °C for 24 hours.

Remove the culture medium.

Add 100 μ L of VSV in MEM.

Incubate at 37 °C for 24 hours.

Remove the culture medium.

Add 50 μ L of 0.05% crystal violet, 10% formalin, 10% ethanol.

Remove the culture medium.

Wash with tap water.

Dry the plate.

Score the CPE.

*At a concentration of 3.6 ng/ml, this antibody neutralizes 10 U/ml of mouse IFN- β to 1 U/ml.

Storage: The solution should be stored below -20 °C.

Repeated freezing and thawing should be avoided.

No preservative substance is used in this product.

Note: Research use only.

Reference: Kawade, Y. J. Interferon Res. 1, 61-70 (1980)



COSMO BIO CO., LTD.

[JAPAN]
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU, TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

[Outside Japan]
2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
e-mail: support@cosmobiousa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobiousa.com

Feb. 26, 2026