



Anti Obestatin (Human) Serum

Cat. No. Y461

Lot No. 1717200424

Description: This antiserum was raised in a rabbit by immunization with a keyhole limpet hemocyanin (KLH) conjugate of synthetic obestatin (human) peptide. The product vial contains 50 µL of the titled antiserum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5 mL) solution. It can be used for immunohistochemistry.

Immunogen: Synthetic Obestatin (human)-KLH conjugate

Host: Rabbit

Amino Acid Sequence of Obestatin (human) ¹⁾

FNAPFDVGIK LSGVQYQQHS QAL-NH₂

Product Form: Lyophilized unpurified serum

Size: 50 µL

Reconstitution: Reconstitute the product with 0.5 mL of 0.01M PBS (pH 7.0) to yield a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN₃ 0.1%).

Storage: The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing- thawing should be avoided.

Suggested Working Dilution Range: Dilution range is 1:5000 for immunohistochemistry. Optimal dilution should be determined by each laboratory for each application.

Specificity (based on non-competitive EIA): Obestatin (human) 100%

Positive Control (Immunohistochemistry): Human stomach

Species Tested (Immunohistochemistry): Human

REFERENCES:

1. Zhang, JV et al: Obestatin, a peptide encoded by the ghrelin gene, opposes ghrelin's effects on food intake. *Science*, **310**, 996-999, 2005



COSMO BIO CO., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <http://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
e-mail: support@cosmobiousa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobiousa.com